English

### 2019-nCoV Antigen Rapid Test Kit (Colloidal gold Assay)

Package specification: 1 test kit, 5 tests kit or 25 tests kit

#### Intended use

It is for self-testing and in vitro diagnostic use only. 2019-nCoV Antigen Rapid Test Kit (Colloidal gold Assay) is an in vitro immunoassay. The assay is for the direct and qualitative detection of SARS-CoV-2 antigen in samples from the human anterior nasal cavity area. This reagent is for home use, testing self-collected nasal swab samples in individuals of 18-65 years old. Anyone over the age of 18-65 should collect sample under the guidance of an adult. Those who are unable to conduct the test on their own should seek support.

#### Principle

In the assay, the immuno-colloid gold technique is used to detect antigens (N-Protein) of SARS-CoV-2. The reagent binding pad is coated with anti-SARS-CoV-2 monoclonal antibodies, each labeled with a colloidal gold marker. A nitrocellulose membrane in the test area of one strip is coated with anti-SARS-CoV-2 antibody. The quality control area within the nitrocellulose membrane is coated with goat anti-mouse IgG antibodies. When tested, the antibodies to SARS-CoV-2 form immune complexes with the N-protein of the virus in the sample being tested. As a result of chromatography, the immune complexes move along the membrane and are captured by the anti-SARS-CoV-2 antibodies coated in the test area. In doing so, a visible red line (L-Line). The free colloidal gold markers or immune complexes move forward and bind specifically to the goat anti-mouse antibody coated in the quality control area. In doing so, they form a visible line (C-Line). If the sample does not contain SARS-CoV-2 antigen, test line will not appear, and only the quality control line (C-Line) appears.

#### Structure of the kit Contents

Test device: There are three different package specifications of 1 or 5 or 25 test cassettes containing immobilized anti-SARS-CoV-2 antibodies labeled with colloidal gold, monoclonal anti-SARS-CoV-2 antibodies coated on 1 line and goat anti-mouse IgG antibodies coated on C line as control. Disposable Sampling Swab: 1 or 5 or 25 pieces. Prepacked Extraction Buffer in Sample Dilution Tube: 500 µl per tube, 1 or 5 or 25 tube. Note: The sample extraction buffer must not be used with a mixed batch.

#### Storage and stability

- 1. Store in dry places and protected from light at 2-30°C. The validity period is 18 months.
2. In general, the kit should be used within 30 minutes after opening the aluminum foil pouch.
3. If the temperature is higher than 30°C or the humidity of the environment is higher than 70%, use the kit as soon as possible after opening the aluminum foil bag.
4. The date of manufacture and the expiry date are printed on the outside of the packaging.

#### Precautions

- Read these instructions before use. The instructions should be read carefully and followed.
• Do not use the test kit or its components after the expiration date.
• The device contains material of animal origin and should be treated as a potential biohazard. Do not use if the bag is open or damaged.
• The test devices are sealed in moisture-proof foil pouches. Inspect each foil pouch before opening. Do not use the test if the foil bag is with holes or damaged or not being completely sealed after the reagent is checked in its proper storage of test reagents or components may result in erroneous results.
• Do not use the extraction buffer if it is discolored or cloudy. Discoloration or turbidity may be a sign of microbial contamination.
• The use should not take any decision of medical relevance without first consulting his or her medical practitioner.
• Do not reuse the test kit.
• After you have done the test, put all parts of the kit in the waste bag. Discard the waste bag in accordance with local regulations.
• All patient specimens must be handled and disposed of as if they were biohazardous. Each specimen must be thoroughly mixed with the Specimen Extraction Buffer prior to testing to ensure a representative specimen prior to testing.
• If specimens and reagents are not brought to room temperature prior to testing, assay sensitivity may decrease. If specimens are collected, stored, and transported incorrectly or improperly, false negative test results may occur.
• Avoid skin contact with the buffer.
• If infection with SARS-CoV-2 is suspected based on current clinical and epidemiological testing criteria recommended by public health authorities, samples shall be collected using appropriate precautions to prevent infection and sent to state or local public health authorities for testing.
• Except in a BSL3 laboratory operating in accordance with BSL3 practice, the virus isolation in cell culture and initial characterization of viral pathogens obtained in cultures of SARS-CoV-2 samples are NOT recommended.

### Quality Control


2019-nCoV Antigen Rapid Test Kit has built-in (procedural) controls. Each test device has an internal standard zone to ensure proper sample flow. The user must verify that the colored line in the "C" region is visible before reading the result. Best laboratory practice recommends testing positive and negative external controls to ensure that the test reagents are working and the test is being performed correctly.

#### Test procedure

- 1. Preparation
a) Adjust the kit to room temperature and make sure the room temperature is 2-30 °C.
b) Open the packaging bag, take out the cassette and put it out in a dry surface.
2. Test procedure
2.1 Place the test set horizontally on the table.
2.2 Add sample.
Place a clean dropper tip on the sample tube. Invert the test tube so that it is perpendicular to the sample hole (S). ADD 3-4 drops to the sample hole (S). Set the timer to 15 minutes.
2.3 Read the result.
Positive samples can be detected in 15 minutes after sample addition. Result is invalid after 20 minutes.

#### Interpretation of the results

NOTE:
1. The color intensity in the test region (T) may vary depending on the concentration of analytes present in the sample. Therefore, any staining in the test region should be considered positive. It should be noted that this is only a qualitative test and cannot determine the concentration of analytes in the sample.
2. Insufficient sample, incorrect procedures, or expired tests are the most likely reasons for control line failure.

#### What should I do after test

If the test result is positive: There is currently a suspicion of a COVID-19 infection. Contact your doctor/general practitioner or the local health department immediately. Comply with local guidelines for self-isolation. To have a PCR confirmatory test performed.
If the test result is negative: Continue to comply with all applicable rules regarding contact with others and protective measures. An infection may also be present if the test is negative. If it is suspected, repeat the test after 1-2 days, as the coronavirus cannot be accurately detected in all phases of an infection.
If the test result is invalid: Possibly caused by incorrect test execution. Repeat the test. If the test results remain invalid, contact a doctor or a Covid-19 test center.

Note: Do not take any decision of medical relevance without first consulting your medical practitioner.

#### Restrictions

- 1 This product is for testing samples of individual nasal swab only.
2 A negative result does not exclude the possibility of COVID-19 infection.
3 The test results of this product are for clinical reference only and should not be used as the sole basis for clinical diagnosis and treatment. Clinical management of patients should be determined in combination with their symptoms and signs, medical history, other laboratory tests (especially pathogen detection), response to treatment, epidemiology, and other information.

#### Sample collection and storage

- 1 Sampling:
Nasal swab: The nasal cavity must be moist. Carefully insert the swab into one nostril of the patient. Insert swab tip 2-4 cm until resistance is encountered. Rotate swab 5 times along nasal mucosa to ensure that both mucus and cells are picked up. Repeat procedure with the same swab in the other nostril to ensure adequate sample is collected from both nasal cavities. Pull the swab out of the nasal cavity.
2 Dip the head of the cotton swab into the extraction buffer tube after sampling. Mix well and squeeze the cotton swab 10-15 times by pressing the tube wall against the cotton swab. Allow the swab to stand for 1 minute to ensure as many samples as possible in the sample extraction buffer. Discard the cotton swab.

- 3 2019-nCoV Antigen Rapid Test Kit can be used to detect both viable and non-viable SARS-CoV-2 viruses.
4 Failure to follow the TEST PROCEDURE and RESULT INTERPRETATION may affect test performance and/or invalidate the test result.
5 Results obtained with this assay, especially in the case of weak test lines that are difficult to interpret, should be used in conjunction with other clinical information available to the clinician.

#### Features

CLINICAL EVALUATION: A clinical evaluation was performed to compare results obtained with the 2019-nCoV Antigen Test Kit and RT-PCR. The comparative PCR was performed by pharyngeal cavity (nasopharynx or oropharynx). The results are summarized below.

Table 1: Antigen rapid test compared to RT-PCR. Clinical diagnosis (PCR results) table with columns for Test results, Positive (+), Negative (-), and Total. Results: Positive (+) 217, 0, 217; Negative (-) 8, 459, 467; Total 225, 459, 684.

2 All samples shall be handled as if they were capable of transmitting infectious agents.
3 Swabs should be tested as soon as possible after collection. Fresh specimens should be used for optimal testing.
4 Do not use specimens that are obviously contaminated with blood, as this may interfere with specimen flow when interpreting test results.

Clinical Sensitivity: 96.44% (93.11%-98.45%)\*
Clinical Specificity: 100% (99.20%-100%)\*
Overall Coincidence Rate: 96.62% (94.98%-97.84%)\*
95% confidence interval

### DETECTION LIMIT

The LoD of Rapid Antigen Test is Ct 30.1.

### PRECISION

Repeatability: Three different QC samples including negative control, weak positive control and moderate positive control were tested using 3 different batches of the product and each sample was tested 6 times obtained on the same testing day by different operators to validate the repeatability of the product.
Reproducibility: Three different QC samples including negative control, weak positive control and moderate positive control were tested using 3 different batches of the product and each sample was tested 10 times obtained on the different testing days by different operators to validate the reproducibility of the product.
The Rapid Antigen Test showed good repeatability on test results by different operators and on the same testing day.

Repeatability: Three different QC samples including negative control, weak positive control and moderate positive control were tested using 3 different batches of the product and each sample was tested 10 times obtained on the different testing days by different operators to validate the reproducibility of the product.
The Rapid Antigen Test showed good reproducibility on test results by different operators and on the different testing days.

### CROSS-REACTIVITY

The Rapid Antigen Test has no cross-reactivities with microbes tested below. Besides, positive test results are not interfered by these microbes.

Table 2: Cross-reactivity. Table with columns: Microbes, Concentration. Lists various microbes like HCoV-229E, HCoV-OC43, HCoV-HKU1, HCoV-NL63, MERS-CoV, Adenovirus, Human Metapneumovirus, etc. with their respective concentrations.

### INTERFERING SUBSTANCES

The following substances, which may be occurring in breath specimens or artificially introduced into the respiratory tract, were evaluated at the concentrations listed below. None of them was found to interfere with the test performance of the 2019-nCoV Antigen Rapid Test Kit.

Table 3: Substances that did not affect the kit. Table with columns: Group, Substance, Final Conc. Lists substances like Mucin, Human blood, Biotin, BEGGI, MOREBENDO, S-O-S, Rhinocort, Rhinasynt, NASONEX, Flixonase, NeilMed, Rhinase, Ayr, Mucinex, HALL'S Relief, LIM, Ribavirin, Oseltamivir, Zanamivir, Amantadine, a-interferon.

Hook effect: Within the titer range of clinically strong positive samples of SARS-CoV-2 antigens, there is no hook effect in the test results of this product.

Ordering information: Beijing Diagreat Biotechnologies Co., Ltd Room 101 Floor 1-3, Building 19, No. 2 Yard, Huanke Middle Road, Tongzhou District, Beijing, China. Website: www.diagreat.com. Email: info@diagreat.com

MedNet EC-REP GmbH Borkstrasse 10, 48163 Munster, Germany
Disposable Sampling Swab manufacturer information: JINAN BABIO BIOTECHNOLOGY CO., LTD. 303, Building 5 of SME Industrialization Base of Huankou Middle Road, Tongzhou District, Beijing, China. Park, 1777 Dazheng Road, High Tech Zone, Jinan City, 250101 Shandong Province, China

INTERTEK SEMKO AB Torshammsgatan 43 Box 1103 SE 164 22 KISTA

Legend for the symbols. Table with columns: Symbol, Description, EU conformity. Includes symbols for Manufacturer, Authorized representative, In vitro diagnostics, Production date, Caution, Protect from moisture, Temperature limited (2-30 °C), and Follow the instructions for use.

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